

# Trillium Brewing: Quality Control Lab

## Canton, Massachusetts




### Featured Products:



CABINETS



EPOXY TOPS



PEGBOARDS



MIXING FAUCETS

Trillium Brewing Company, a farmhouse-style craft brewery, needed a quick fix for their small quality control laboratory at their second location in Canton, Massachusetts. New England Lab partnered with the Rapid Ship Division, NEL Direct, to supply and install casework and furniture within four weeks. With this small room setback from the merchandise and tasting bar, the installation proceeded without disruption during retail operation hours. NEL Direct's Bright White painted steel casework, and adjustable height tables were provided for the project based on the client's requirements. Black epoxy resin work surfaces, a phenolic pegboard, plumbing fixtures, and sinks completed the space.

